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c.0 m.0 y.0 k.80
Pantone 425 C

Material

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Papel: Sulfite 63g/m²

Dobra: Vinco vertical

Tamanho: 90 x 140mm

Marcas de Corte/Dobra

———— Corte Linear

- - - - - Vinco de dobra

Resinlay Pattern Photo Gel

(PT)

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Light cured universal modelling acrylic resin in gel.

Resinlay Pattern Photo is a light-curing acrylic resin developed for universal modelling and casting patterns. Available in gel and paste consistencies, both versions can be used in associated techniques, allowing a wide number of applications by a complete polymerization and burning clearly leaving no residues, even on thicker layers. Its gel consistency allows an easy and precise application, using the application tips, due to its thixotropic characteristic.

Indication: This product is indicated for use by both, the dentist and the dental laboratory technician. Among its main indications are the manufacture of: **Dentist:** Fixation and union of transfers in cases of multiple prostheses on two or more implants; Customization of transfers for copying the emergency profile; Copings pattern sculpture; Copings for molding; Fixation of temporary splint for soldering on bars for protocols or fixed prosthesis structures; Intra canal impression for post and Core by direct technique; Fabrication and adjustment of guides: surgical, for positioning and insertion of provisional crowns; Reliefs. **Prosthesis Technician:** Preparation of body test for: implant bars, veneers, prostheses structures, copings; Abutments; Post and Core fabrication for casting, which can be used in association with paste version; Caps for molding; Reliefs.

Use: Prior to application, the area shall be isolated with the product Isolar TDV or glycerin. Remove the syringe cap and attach the applicator tip. Apply the resin directly on the area to be worked on. It is suggested to the application be performed by layers of up to 2mm thickness. In case of application of a very thick layer, it is recommended to apply thinner overlapping intermediate layers. Before building the next layer, light cure the surface for 20 seconds, without removing the dispersion layer. After the use, the cap should be replaced in the syringe, preventing the resin from polymerizing. It is recommended to remove the residual material from the applicator tip with a jet of compressed air, ensuring that the resin does not polymerize and obstructs the resin flow in further applications. Perform the finishing of the piece. The dispersion layer can be removed by the finishing. For intra-oral application, the material is for

transitory use, that is, it can be applied in the mouth, only during the molding and sculpture process.

Warnings, precautions and restrictions: Professional use. The product must not be used in patients who show allergic reactions to any of its components. To prevent the product from polymerizing, close the syringe immediately after use. The packaging should be protected from heat and direct exposure to sunlight. The Material Safety Data Sheet for this product can be obtained by telephone. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Cleaning: Before sterilization, execute the cleaning with the aid of water and enzymatic detergent.

Sterilization: The applicator tips are autoclavable (126°C, 16min) after each patient. The sterilization can cause slight change in the coloring, but does not affect its functionality.

Storage: To be stored in a cool dry place, protected from sunlight. Do not store above 30°C.

Composition: Methacrylate and activators.

Disposal: Waste must be disposed of in an industrial or sanitary landfill, in compliance with local legislation in force. Do not dispose of waste in sewers, waterways, or household waste bin.

| Ref. | Código/Code | Conteúdo/Contenido/Content |
|--------|-------------|--|
| 0.1260 | 0012600002 | 3g de gel + 3 ponteiras / 3g de gel + 3 puntas / 3g gel + 3 tips |

Manufacturer: TDV Dental Ltda.
CNPJ 81.591.786/0001-60 | XV de Novembro, 9944
CEP 89107-000 | Pomerode/SC - Brasil
Reg. ANVISA 10291220100
☎ 55 47 3395-6115 | 0800 047 1020
info@tdv.com.br | www.tdv.com.br

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